

New from
Katena



katena products, inc.
4 Stewart Court
Denville, NJ 07834
☎ 973-989-1600
☎ 973-989-8175
e-mail globe@katena.com
www.katena.com



Chang "Quick Chop" Phaco *Combination Instrument*



K3-2368 Chang Combination Chopper, for right hand use
K3-2369 Chang Combination Chopper, for left hand use

This double-end instrument features a modified Lieberman Microfinger and a pointed Maloney "Quick Chop" tip. The Microfinger tip, used for horizontal chopping of softer nuclei, is more delicate than the original with a defined inferior edge for efficient dividing. The "Quick Chop" tip, used for denser nuclei, is conical with a sharp point to easily penetrate and split the nucleus vertically.

Developed with David Chang, MD of Los Altos, CA